

Methoxy polyethylene glycol-epoetin

beta (Mircera[®]) Protocol

Methoxy polyethylene glycol-epoetin beta (Mircera®) ICD 10 code D63.1

Anemia in chronic kidney disease

Purpose: To provide optimal management of ESRD related anemia in dialysis patients

Hemoglobin Target Goal: 10.0-11.0 g/dL

Methoxy polyethylene glycol-epoetin beta Dosing:

Doses are based on estimated dry weight and rounded to the following steps:

Step	Dose
1	30 mcg every <i>four</i> weeks
2	50 mcg every <i>four</i> weeks
3	30 mcg every two weeks
4	50 mcg every two weeks
5	60 mcg every two weeks (30 mcg + 30 mcg)
6	75 mcg every two weeks
7	100 mcg every two weeks
8	150 mcg every two weeks
9	200 mcg every two weeks

Table 1

- 1. Methoxy polyethylene glycol-epoetin (Mircera®) will be increased and decreased in 1-step or 2-step increments, based on scale above.
- 2. Mircera® will be administered IV to in-center hemodialysis patients, and SQ to home dialysis patients.
- 3. Mircera® ceiling is 200 mcg every two weeks (or 3.0 mcg/kg/2 weeks, whichever is lower). Orders above 200mcg every two weeks require facility medical director or CMO approval.

Initiating Mircera® for new patients or ESA naïve patients

For new patients or established patients who have not received an ESA within the last 3 months, initiate as follows:

- 1. Iron repletion per iron standing orders
- 2. AND
 - a. If Hgb < 10 g/dL, then start Mircera[®] at 0.6 mcg/kg/2 weeks, and round down to closest step per Table 1 but no less than 30 mcg every 2 weeks (Step 3).
 - b. If Hgb 10.0-10.4 g/dL, then start Mircera® at 30 mcg every 2 weeks (Step 3).
 - a. If Hgb >= 10.5 g/dL, then do not start Mircera® until patient meets criteria.

Patient Name _	NKC#
MEC approved 9.12.19	Page 1 of 3

Mircera® Dosing Adjustment

1. Titrate Mircera[®] per the following table for patients who have a Mircera[®] order and had not been changed in the last 4 weeks:

Mircera® Dosing Adjustment					
Hgb decreased by greater than or equal to 0.5 g/dL since last dose change					
Hgb (g/dL)	Step Dose Change				
Less than 10	2 step dose increase				
10.0-10.9	1 step dose increase				
11-11.9	No Change				
Hgb increased/decreased by less than 0.5 g/dL since last dose change					
Hgb (g/dL)	Step Dose Change				
Less than 9.5	2 step dose increase				
9.5-9.9	1 step dose increase				
10.0-10.4	1 step dose increase, if Hgb decreased; do not change if Hgb				
	increased or stayed the same				
10.5-11.4	No change				
11.5-11.9	1 step dose decrease; if patient is on Step 1, do not HOLD				
Hgb increased greater than or equal to 0.5 g/dL since last dose change					
Hgb (g/dL)	Step Dose Change				
Less than 10	1 step dose increase				
10-10.4	No Change				
10.5-11.9	1 step decrease; if patient is on Step 1, do not HOLD				
Hgb (g/dL)	Dose Change				
Greater than or	Hold Mircera; check Hgb at next redraw for home dialysis				
equal to 12 g/dL	patients, and every week for in-center patients.				
If Hab is increased or decreased at least 1.0 a/dl since the last Hab level:					

If Hgb is increased or decreased at least 1.0 g/dl since the last Hgb level; recheck Hgb within next 2 dialysis treatments for in-center HD and at next redraw for home patients. Follow the algorithm based on the results of the recheck, e.g., if the value remains the same as the first draw, then follow the algorithm for no change.

Table 2

- 2. Do not change Mircera® dose more frequently than every 4 weeks EXCEPT:
 - a. If Hgb falls from above 10 g/dL to less than 10 g/dL, increase dose after 2 weeks.
 - b. If Hgb is already less than 10 g/dL and drops greater than 0.5 g/dL, increase dose after 2 weeks.
 - c. If Hgb >= 12 g/dL, hold Mircera[®] and check Hgb every week for in-center patients, and at next redraw for home dialysis patients. Resume Mircera[®] with 1-step decrease as soon as Hgb is < 11.8 g/dL and last dose was administered 2 weeks ago or more. If Hgb remains >= 12 g/dL for more than 2 months, return to regular Hgb testing policy.

Conversion from darbepoetin or erythropoietin to Mircera®

- 1. When a patient with a darbepoetin (Aranesp) or erythropoietin order switches to Mircera[®], discontinue darbepoetin (Aranesp) or erythropoietin order.
- 2. Convert darbepoetin or erythropoietin to appropriate dose of Mircera[®], per conversion dose chart below. Convert to Mircera[®] when the next ESA dose is due.

MEC approved 9.12.2019 Page 2 of 3

3. If ESA is on HOLD from another protocol, wait until Hgb is less than 11.8g/dl, then convert ESA as follows: See Table 3 or 4 to convert previous ESA dosing to Mircera® Step, then see Table 1 and decrease 1 Step.

Erythropoietin to Methoxy Polyethylene Glycol Epoetin- beta Conversion Dose Chart				
Epogen Dose (U) per week - total	Mircera [®] Dose			
	Dose (mcg)	Frequency		
< 2000	30	Every 4 weeks		
2000 - < 3000	50	Every 4 weeks		
3000 - < 5000	30	Every 2 weeks		
5000 - < 8000	50	Every 2 weeks		
8000 - < 11,000	60	Every 2 weeks		
11,000 - < 18,000	75	Every 2 weeks		
18,000 - < 27,000	100	Every 2 weeks		
27,000 - < 42,000	150	Every 2 weeks		
>= 42,000	200	Every 2 weeks		

Table 3

Darbepoetin (Aranesp) to Methoxy Polyethylene Glycol Epoetin- beta Conversion Dose Chart				
Darbepoetin Dose (mcg) per week - total	Mircera [®] Dose			
	Dose (mcg)	Frequency		
< 10	50	Every 4 weeks		
10 - <20	30	Every 2 weeks		
20 - <30	50	Every 2 weeks		
30 - < 40	60	Every 2 weeks		
40 - < 50	75	Every 2 weeks		
50 - < 60	100	Every 2 weeks		
60 - < 100	150	Every 2 weeks		
>= 100	200	Every 2 weeks		

Table 4

Labs:

Draw CBC per routine lab orders.

Suzanne Watnick, MD Physician Name (Please Print)	
Physician signature	Date
Patient Name	NKC#