

# Chronic Maintenance In-Center Hemodialysis Standing Orders (Addendum)

---

Chronic Maintenance In-Center Hemodialysis Standing Orders – Medications PRN,  
item 11b clarification:

11b.

## Staff Application

Apply a moderate layer of lidocaine 2.5% / prilocaine 2.5% cream (or other generic equivalent) at the site of needle placement 15-60 minutes (preferably 30 or greater minutes) prior to start of hemodialysis.

## Patient Application

Dispense to patient: One 30-gram tube of lidocaine 2.5%/prilocaine 2.5%.

Directions for use: Apply a moderate layer (approximately 2.5 grams) of cream at site of needle placement 15-60 minutes (preferably 30 minutes or greater) prior to start of hemodialysis. Cover area with occlusive dressing (plastic wrap).

Refills: As needed. Not to exceed one 30-gram tube per 12 dialysis treatments.

Matthew Rivara MD  
Physician Name (Please Print)

\_\_\_\_\_  
Physician signature  
(*see Initial Orders*)

\_\_\_\_\_  
Date

**Patient Name:** \_\_\_\_\_ **NKC#** \_\_\_\_\_